


ENDFLO™ by TechDevice

The ENDFLO™ catheter is a single use device designed for temporary vessel occlusion. Select from balloon size and shaft working length options below. The product is made from Pebax and Latex.



Part Number	Balloon Size (mm)	Working Length (cm)	Shaft Size
Endflo 5001	27	45	5 French
Endflo 5002	27	65	5 French
Endflo 5003	27	100	5 French
Endflo 5004	18	45	5 French
Endflo 5005	18	65	5 French
Endflo 5006	18	100	5 French
Endflo 5007	15	45	5 French
Endflo 5008	15	65	5 French
Endflo 5009	15	100	5 French
Endflo 5010	11	45	5 French
Endflo 5011	11	65	5 French
Endflo 5012	11	100	5 French

- Configuration: Two lumens, one (1) for inflation/deflation, 2nd for guidewire
- Indications: For temporary vessel occlusion where stoppage or blood is desired in areas other than neuro and cardiac use
- Approvals: TechDevice 510(k) # K051137 for temporary vessel occlusion
- Packaging options include: Bulk non-sterile, individually packaged in trays, pouched, labeled and sterilized

6/15/10

TechDevice
CORPORATION

ISO 13485
CERTIFIED

1-888-TECHDEV • www.techdevice.com • 617-972-5800

FDA
REGISTERED