

TDXPress HEVEA LATEX BALLOONS

Shorten development cycles with TechDevice Hevea Latex Balloons.

Minimum order of 10 units.

- 2.25" lengths
- wall thickness of .005" - .008" or .008" - .011"
- inside diameters made to fit standard French sizes
- one round closed end



Items in stock can ship either same or next business day. If an item is not on the shelf, your sales person can provide you with a lead time for the quantity you need.

TechDevice manufactures custom balloon sizes, shapes and lengths. Ask your sales person for details or provide a specification for us to review.

Latex Sleeves Shaft Size	Part Number	
	.005" - .008" Wall Thickness	.008" - .011" Wall Thickness
.037"/3 FR	40118-01	40118-14
.053"/4 FR	40118-02	40118-15
.066"/5 FR	40118-03	40118-16
.079"/6 FR	40118-04	40118-17
.092"/7 FR	40118-05	40118-18
.105"/8 FR	40118-06	40118-19
.118"/9 FR	40118-07	40118-20
.131"/10 FR	40118-08	40118-21
.144"/11 FR	40118-09	40118-22
.158"/12 FR	40118-10	40118-23
.170"/13 FR	40118-11	40118-24
.184"/14 FR	40118-12	40118-25
.197"/15 FR	40118-13	40118-26

- All balloons are shipped bulk non-sterile
- Lead time may vary – check at time of order

12/1/11

TechDevice Stock Products – Short Lead Times
CORPORATION



1-888-TECHDEV • www.techdevice.com • 617-972-5800

