


# TDXPress URETHANE BALLOONS

| Part Number | Material | Dia. (A) (mm) | Balloon Length (B) (mm) | Overall Length (mm) | Prox/Distal ID (inches) | Prox/Distal Cones (degrees) | Rated/Avg. Burst (atm) |
|-------------|----------|---------------|-------------------------|---------------------|-------------------------|-----------------------------|------------------------|
| 40131-02    | Urethane | 2.0           | 8                       | 13.5                | .035/.035               | 20/20                       | 12/16                  |
| 40132-12    | Urethane | 2.0           | 20                      | 25.5                | .035/.035               | 20/20                       | 12/16                  |
| 40132-13    | Urethane | 2.0           | 30                      | 35.5                | .035/.035               | 20/20                       | 12/16                  |
| 40131-06    | Urethane | 3.0           | 15                      | 23.25               | .039/.039               | 20/20                       | 12/16                  |
| 40131-14    | Urethane | 3.0           | 20                      | 28.25               | .039/.039               | 20/20                       | 12/16                  |
| 40131-15    | Urethane | 3.0           | 40                      | 48.25               | .039/.039               | 20/20                       | 12/16                  |
| 40127-09    | Urethane | 4.0           | 20                      | 35.94               | .093/.093               | 34/34                       | 10/14                  |
| 40127-08    | Urethane | 4.0           | 30                      | 25.94               | .093/.093               | 34/34                       | 10/14                  |
| 40127-01    | Urethane | 4.0           | 40                      | 45.94               | .093/.093               | 34/34                       | 10/14                  |



Urethane is a semi-compliant balloon material that:

- is used for lower pressure applications
- is typically sterilized by means of Eto
- can be thermally welded or adhesive bonded

6/15/10

**TechDevice** Standard Products – Short Lead Times  
CORPORATION



ISO 13485  
CERTIFIED

617-972-5800 • www.techdevice.com

FDA  
REGISTERED